



R82[®]


Flamingo

Adjustable, functional and reuseable

Why Flamingo?

Easily adapted, functionally superior, cost effective to reuse, offering both low and high level support.

Features



Which Flamingo?

- What do the child need to do when in the seat?
- Toileting? bathing? Combination of both?
For how long?
- What support is needed in the seat?
- Can the child sit Independently? Is more support needed at the hips, trunk, head?
- Where is the seat going to be used?
- Free standing in the bath, over the toilet, in the shower or a combination of these?



Choosing a base

Our Flamingo family is one of the most versatile products on the market with up to five configurations available. It can be used freestanding with a pan, directly over the toilet, as a shower chair and as a potty chair. The seat is also removable for direct mounting on a toilet or in a bath.



Suction cups fitting to bath (R82 Orca)



Potty frame adaption



Direct fitting to the toilet seat



Use with a commode pan



Flamingo frame



Flamingo high-low frame

From a therapist perspective

Flamingo provides possibilities for staged development of toileting and hygiene skills. Helena Poulton, OT & Independent Consultant gives her opinion:

The right Flamingo for the child

- When using as a toilet seat, support can be fully adjusted to find the optimum position of toileting by changing the angle and position of the hips, knees and feet
- Support through a contoured seat and 4 point hip belt. Can be integrated with accessories
- Supportive, adjustable support to promote alignment and comfort so the child is safe and can enjoy bath/shower time
- The tilt and height adjustment allows for much safer and easier hoisting

The right Flamingo for carers

- Reduces number of transfers needed between equipment. A total hygiene solution in one chair. Easy to clean
- Easy to adjust while the child is in the seat between bathroom tasks where the hip angle and tilt in space often need changing for function
- Responsive foot pedals to adjust height & tilt in space
- Excellent manoeuvrability even in tight spaces

From a parents perspective

The Flamingo High-Low's height adjustability function allows Tia's parents to perform bathing and hygiene tasks at a convenient and good working position.

Tia's mum, Lesley tells us:

"Not only does the chair help me but it will aid the development of Tia's independent bathing and hygiene skills...We have a smooth transition from shower to toilet to sink"



New and improved Flamingo High-low

One of the most adaptable product for toileting & bathing in the market. Flamingo High-low is a modern and comfortable toilet/bathing chair that fits into every bathroom.

Features

Armrests/Tray

Padded tray offers arms support and swings away for transfers

Head support

Height and angle adjustable with 4 different mould shapes to suit individual need

Splash guard

Soft removeable splash guard, also helps as a guide for midline positioning

Side supports


Fixed or Swing away supports. Swing away option is provided to make washing easier when showering. Curved supports - shaped for maximum comfort and support

Footsupport

Adjustable for optimum positioning and comfort

Compatibility

The new frame design means Flamingo will go lower to the toilet and fit an increased range of modern toilets



Read more about the
R82 Flamingo High-low on our website...

Accessories


Numerous accessory options can be chosen for the Flamingo seat. You can explore some of the accessories below and see the complete Flamingo seat and its range of accessory options on our website.



Back extension



Head supports




Side supports



Armrests/tray



Sides with upholstery



Splash guards



Commode pan



Vests and belts

Early Intervention

The Penguin bath seat provides head and trunk support and freedom of movement for arms and legs. It has an adjustable abduction block to fit the individual child.

Flamingo seat, size 1 & 2 fits in the Orca bath for a temporary solution for the safe bathing of young children. Frame can be removed for easier storage.



Other related products

Mobile pediatric toilet and bath chairs shall provide children, and youngsters with special needs the best support and the carers a flexible, optimal working conditions in any environment.

R82 Swan Curo



The Swan Curo is a mobile toilet and bath chair, ensuring everyday convenience for user and caregiver.

The shell seat is carefully designed to enclose the user for high comfort and has smooth surfaces, making it easy for the caregiver to clean. The shape of the shell seat along with the soft PU inlay, ensures high seating comfort for the user. Positioning of the user is easily done using the tilt function assisted by a gas spring and the height- and angle-adjustable foot support. To obtain the best possible height for transfer, for caregivers or the toilet, the frame is height-adjustable.

R82 Manatee



The Manatee is constructed to provide both caregivers and users a wonderful bathing experience.

The bathing chair is designed to support the user to maintain safe sitting position. The depth of the seat can be adjusted, allowing the bathing chair to “grow” with the child and provides many years of use. The adjustments and accessories are logically positioned, affording the carer better working conditions. The Manatee is very lightweight and can easily be folded for storage or transportation.

R82 Heron



The Heron is a state of the art toilet and bathing seat.

The simple design and clean appearance feature multiple functionalities that improve both the user and carer experience. It sets new standards for users with special needs for positioning. The seat is ergonomically correct, meaning it can be tilted and the back reclined without any disadvantage for the user. Clever design means that during recline, the anatomical pivot point at the axis of the hips is maintained. This provides the user more comfort when seated.



Flamingo_Brochure B1600_01_INT | Reserving our rights for printer's error or discontinued products.
Please find an updated edition on our website etac.com © R82 A/S All rights reserved.

Etac is a world-leading developer of innovative assistive devices and patient handling equipment. Since 1973, we have been committed to improving quality of life for the individual, family members and caregivers.

For the intended purpose of the products and latest information, please visit www.etac.com

R82 A/S

+45 79 68 58 88

R82@etac.com

www.etac.com

etac[®]
Creating Possibilities