



**K300**  
PS JUNIOR

## Ready? Go!

### Playing on wheels

The **Permobil K300 PS Junior** enables children to join in all the fun. Its fully adjustable seat tray allows for the child's growth. The adjustable armrests, electric seat lift and tilt function, programmable R-Net controls and the UniTrack system offer extra comfort. This electric wheelchair will withstand the odd knock or bump and is available in stylish and vibrant colours. Watch the video on our site!



**permobil**



## Lighting

Optionally available with energy efficient front and rear LED lighting.



## Control

R-net electronics: intelligent control system (ICS) with programmable settings.



## Push handles

Push handles come as standard with the K300 PS Junior.



## Seat tilt

45° powered or manual seat tilt adjustment is available as an option.



## Armrests


Easily-adjustable armrests for convenient adjustment of height, angle and width.



## Footrests

Footrests rotate away in one move, with adjustable and extremely strong footplates with several axes. One-piece leg support also available.


## Colours



cobalt blue




black onyx



phantom gray



popstar pink



deep purple




volcano red

## Upholstery



leatherette



black fabric

## Specifications

Total length	1080 mm
Total width	620 mm
Minimum transport length	885 mm
Seat height (fixed)	460 mm
Seat height with powered seat lift	460 - 660 mm
Total weight incl. batteries	109 kg
Maximum user weight	76 kg
Seat width	275 - 400 mm
Seat depth	250 - 450 mm
Backrest width	275 - 400 mm
Backrest height	350 - 560 mm
Distance between armrests	adjustable
Backrest angle - adjustable	85° - 120°
Leg support angle - adjustable	-

Battery capacity	2 x 60 Ah
Range (ISO 7176-4)	30 - 35 km
Maximum speed	7 km/hr
Turning radius (ISO 7176-5)	625 mm
Obstacle height	60 - 80 mm
Suspension	yes
Electronics	R-net
Powered seat lift	yes
Powered seat tilt	yes
Powered leg rest adjustment	no
Powered backrest adjustment	no
Guarantee	1 - 2 years
ISO 7176-19/ISO 10542-3	yes - crash tested
EN12182/EN12184 CE	yes



Τηλ.: 77771314, 24642070, Κινητό: 99455233

[www.medicomstore.com](http://www.medicomstore.com)

**permobil**  
THE POWER OF MOBILITY